# **mott** corporation

## All-Metal Steam Filters for Biotechnology Applications



In most Bioprocess applications, vessels, piping and filter housings are sterilized by steam-in-place (SIP) operations. Many applications also call for steam jacketing of vessels, filter housings and process lines to maintain process temperature.

The steam required to sterilize or operate these systems often contain a considerable amount of particulate in the form of rust or other pipe-scale. If contaminated by these impurities, it diminishes the life of the valves, filters, steam jackets and other ancillary equipment.

Mott High Purity Steam Filters prevent SIP systems from being contaminated and help increase the life of adjoining equipment.

Mott High Purity Steam Filters incorporate sintered porous stainless steel filter media that ensures the effective retention of particles. The 316L SS filter elements provide excellent chemical resistance and mechanical stability due to their special design. Using Mott High Purity Steam Filters can considerably reduce the post-filtration costs of your system.

## Cartridge Features

- Porous sintered stainless steel construction
- High flow rate with low differential pressure
- High mechanical stability
- High temperature range

### Filter Cartridge Specifications

- Available in media grades 10, 20, and 40 μm
- Welded end caps either flat end or centering cone
- 226 with locking flanges or 222 adapters with double "O" ring

Nominal Dimensions	10 inch Cartridge	20 inch Cartridge	30 inch Cartridge	40 inch Cartridge
Overall Length	10.75 Inch/273 mm	20.75 Inch/527 mm	30.75 Inch/781 mm	40.75 Inch/1035 mm
Outside Diameter	2.5 Inch/63.5 mm	2.5 Inch/63.5 mm	2.5 Inch/63.5 mm	2.5 Inch/63.5 mm
Filtration Area	0.54 Ft <sup>2</sup> /0.05 m <sup>2</sup>	1.09 Ft <sup>2</sup> /0.10 m <sup>2</sup>	1.63 Ft <sup>2</sup> /0.15 m <sup>2</sup>	2.18 Ft <sup>2</sup> /0.20 m <sup>2</sup>
Filter Media	316L SS	316L SS	316L SS	316L SS
End-caps, Adapter	316L SS	316L SS	316L SS	316L SS
Maximum Forward Differential Pressure in Air, Oxygen or Nitrogen	250 psid (17.23 Bar)	250 psid (17.23 Bar)	250 psid (17.23 Bar)	250 psid (17.23 Bar)
Maximum Reverse Differential Pressure in Air, Oxygen or Nitrogen	250 psid (17.23 Bar)	250 psid (17.23 Bar)	250 psid (17.23 Bar)	250 psid (17.23 Bar)
Maximum Operating Temperatures	450°C (842°F) @ 125 PSIG			
Cartridge Connections Double O-Ring	226 with Locking Flange or 222			
End caps	Flat end cap or centering cone end cap			

## high purity products



222 Double O-Ring Cartridge

#### 226 Double O-Ring Cartridge with Locking Flanges



## Bio-Pharm Steam Filter Cartridges Part Number

	Code	Product	Family					
	6840	Steam Fil	ter Elemen	it				
		Code	Connect	or Type				
	ĺ	222	Double C	<b>D-Ring Size</b>	222 Conne	ector		
		226F	Double C	<b>D-Ring Size</b>	226 Conne	ector with L	ocking F	langes
	ĺ		Code	Porous	Tube Leng	th (Dimens	ion A at	pove)
	ĺ	ĺ	10	Porous 7	lube Length	n, 2.5 Inch (	DD	
			20	Porous 7	Tube Length	n, 2.5 Inch (	DD	
			30	Porous 7	lube Length	n, 2.5 Inch (	DD	
		ĺ	40	Porous 7	Tube Length	n, 2.5 Inch (	DD	
		Ì		Code	Media Gr	ade		
				10	10 µm Me	edia Grade	Porous	
				20	20 µm Me	edia Grade	Porous	
				40 	40 µm Me	edia Grade	Porous	
					Code	End Cap		
					A00	316L SS I	Flat End	Сар
					F12	3/4 Inch E	Diameter	, 60° Centering Cone
						Code	Porous	s Media
						A	316L S	S Porous Media
							Code	Hardware Material
							А	316L SS Hardware
Example:	6840	226F-	10-	0.2-	A00-	Α	Α	

### Order O-Ring Separately

Code	O-Rings			
2-222	Standard	222 size O-rings		
2-226	Standard	226 size O-rings		
	Code	Materials		
	63	Neoprene		
	64	Teflon®		
Í	65	Viton		
İ	66	Buna-N		
İ	71	Ethylene Propylene (EPDM)		
İ	72	Silicone		
İ	86	Teflon <sup>®</sup> Encapsulated-Silicone		
İ	87	Teflon <sup>®</sup> Encapsulated-EPDM		
İ	91	Kal Rez		
İ	92	Teflon <sup>®</sup> Encapsulated Viton		
İ	1A	Viton Steam Resistant		

Example: 2-222- 63

TEESING WE MAKE YOUR TECHNOLOGY WORK

 $\textit{Teflon}^{\circledast}$  is a registered trademark of E. I. DuPont Nemours & Co., Inc.

## **mott** corporation

